




HAINES
EDUCATIONAL



Vernier[®]
SCIENCE EDUCATION



Haines Educational & Vernier

Your Partners in Datalogging Solutions



Deepen Student Engagement with Our Software

Vernier Graphical Analysis® Do More with Data

Easily visualise and interact with experiment data from nearly any Vernier sensor to help students make sense of the phenomena they can't always see. Essential tools are free - unlock advance features like data sharing and enhanced analysis with Pro.

Learn more at: haines.com.au/vernier/graphical-analysis

Vernier Video Analysis® Study Motion Everywhere

From a basketball's arc to a bottle rocket's launch, this intuitive app makes it easy for students to turn videos into dynamic motion investigations.

Learn more at: haines.com.au/vernier/video-analysis

Vernier Data Collection and Analysis Bundle Bundle and Save

Get a 3-year subscription to Graphical Analysis Pro and Video Analysis, alongside easy access to our free Spectral and Instrumental Analysis apps.

Learn more at: haines.com.au/vernier/dca-3yr

Vernier Connections® Phenomenon-Based, 3D Learning

Explore web-based, interactive lessons designed to put hands-on investigation at the heart of students' learning.

Learn more at: haines.com.au/vernier/connections



LabQuest 3

VERLQ3

LabQuest 3 - a powerful, advanced, easy-to-navigate and versatile data-logging solution for STEM students. The built-in data collection and analysis app works with all Vernier sensors, supporting hands-on data collection in the classroom, in the lab and in the field.




Scan QR Code for more information

Go Direct[®] Temperature Probe

MOST POPULAR

VERGD001

Go Direct Temperature Probe - The Go Direct Temperature Probe is a rugged, general-purpose sensor that students can use to monitor temperature.



Scan QR Code for more information

Go Direct[®] Sensor Cart

VERGD022 & VERGD021

The Go Direct Sensor Cart features built-in, wireless sensors for one-dimensional dynamics and kinematics experiments.



VERGD022

VERGD021




Scan QR Code for more information

Go Direct[®] Colorimeter

VERGD006

The 4-wavelength Go Direct Colorimeter measures absorbance or transmittance of a liquid sample.



Scan QR Code for more information

FOR ALL OTHER VERNIER PRODUCTS, PLEASE VISIT OUR WEBSITE.

Go Direct[®] pH Sensor

VERGD002

The Go Direct pH Sensor is a general-purpose pH sensor used to monitor the pH of aqueous solutions.



Go Direct[®] Force and Acceleration Sensor

VERGD019

The Go Direct Force and Acceleration Sensor couples a 3-axis accelerometer with a stable and accurate force sensor that measures forces as small as ± 0.1 N and up to ± 50 N. Measure pushes and pulls in the classroom or outdoors.



Go Direct[®] CO₂ Gas Sensor

MOST POPULAR

VERGD009


The Go Direct CO₂ Gas Sensor measures gaseous carbon dioxide concentration levels, air temperature, and relative humidity. It can be used in a variety of biology experiments.



Go Direct[®] Motion Detector

VERGD020

The Go Direct Motion Detector uses ultrasound to measure the position, velocity, and acceleration of moving objects.



HOW TO ORDER:














www.haines.com.au



sales@haines.com.au

The best selling Go Direct[®] Sensors by Subject

PRODUCT	ORDER CODE	 BIOLOGY	 CHEMISTRY	 PHYSICS	 ENVIRONMENTAL SCIENCE	K - 8 K-8
Temperature Probe 	VERGD001	●	●	●	●	●
Sensor Cart 	VERGD022 (GREEN) VERGD021 (YELLOW)			●		
Colorimeter 	VERGD006	●			●	●
pH Sensor 	VERGD002	●	●		●	●
Force and Acceleration Sensor 	VERGD019			●		●
CO₂ Gas Sensor 	VERGD009	●	●		●	●
Motion Detector 	VERGD020			●		●

* View the complete Vernier range at www.haines.com.au

Haines Educational & Vernier

Haines Educational proudly stocks the **entire range of Vernier data loggers, sensors, and probes** – everything you need to bring hands-on, data-driven science to your classroom or lab.

Whether you're upgrading your existing setup or starting fresh, our team can help you choose the right Vernier solutions for your curriculum and budget.

Contact us today for expert support, best pricing and access to the complete Vernier collection.



Contact us now for more information:

Haines Educational

1/82 Lewis Road

Wantirna South VIC 3152

Phone: 1300 330 232

Email: sales@haines.com.au

www.haines.com.au



Go Direct®

All Go Direct®
Sensor Accessories
connects via
Bluetooth® Wireless
Technology or
via USB to your
devices.

